FORA 6 Total Cholesterol Test Strip

312-4640100-007 ver 4.0 2021/08

Warnings

- For **In vitro** diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that has come into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your FORA 6 Plus / Connect Multi-functional Monitoring System Owner's Manual before you use this test strip. Use only FORA 6 Test Strips with FORA 6 Plus / Connect Multi-functional Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- Please do not use FORA 6 Plus / Connect Multi-functional Monitoring System on critically ill patients. The collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological total cholesterol level might not be a true reflection.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

FORA 6 test strips, when used together with FORA 6 Plus / Connect Multi-functional Monitoring System is intended for the quantitative measurement of total cholesterol in fresh capillary whole blood from the finger. It is for in vitro diagnostic use only.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use ONLY heparin for anticoagulation of whole blood samples. Please do NOT use EDTA for anticoagulation.

Contents

- Total Cholesterol Test Strip ► Code strip ► Instruction for use

Limitations

- Hematocrit: The hematocrit level is limited to between 20% and 60%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: **This test strip must not be used for the testing of newborns.**
- This test strip is used for testing fresh capillary and venous.
- Altitude Effects: Altitudes up to 2,438 m (8,000 ft) do not affect test results.
- There are no significant interference in the presence of methyl DOPA, acetaminophen, uric acid, dopamine, gentisic acid, ascorbic acid, ibuprofen, salicylate (when at physiological or therapeutic levels).
- Lipemic Effects: Glucose levels up to 476 mg/dL (13.89 mmol/L) do not affect the results significantly.

Storage and Handling

- Use each test strip immediately after taking it out of the vial or individual foil pack. Close the vial immediately after taking out a strip.
- **IMPORTANT: Do not use the test strips if they have expired.**
- Keep the vial closed at all times.
- Store the test strips in their original vial **ONLY**. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Do not bend, cut, or alter the test strip.
- Store the test strips in a cool, dry place between 2°C and 30°C (35.6°F and 86°F) and below 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.

Testing Your Total Cholesterol

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

Your total cholesterol readings deliver plasma equivalent results and is displayed in milligrams of total cholesterol per deciliter of blood (mg/dL).

The measurement range of this meter is 100 mg/dL ~ 400 mg/dL (2.6 mmol/L ~ 10.4 mmol/L).

Reference values:

Total Cholesterol* < 200 mg/dL (5.17 mmol/L)
*Centers for Disease Control and Prevention, CDC (2008)

Please consult your doctor to determine a target range that works best for you.

Questionable or Inconsistent results

- If your test results are unusual or inconsistent with how you are feeling:
- Make sure the confirmation window of the test strip is completely filled with blood.
 - Check the expiry date of the test strips.
 - Check the performance of your meter and test strip with the control solutions.
 - Make sure your monitor has correct coding, and the code is the same as the code printed on the vial label or individual foil packet you are using.

Please Note: If most of your results are unusually high or low, please contact your healthcare professional.

Quality Control Testing

Our control solutions contain a known amount of total cholesterol that can react with test strips. You can check the performance of the meter, test strip and your technique by comparing the control solution results with the range printed on the label or individual foil pack. Checking regularly can ensure your test results are accurate. Please refer to the Owner's Manual for complete testing instructions.

IMPORTANT: The reference range of the control solutions may vary with each new test strip. Make sure you check the range on the label of vial or individual foil pack of your current test strip.

Chemical Components

> Enzyme: Cholesterol esterase (Pseudomonas sp.) ≥ 0.1U
> Enzyme: Cholesterol oxidase (Streptomyces sp.) ≥ 0.1U
> Electron shuttle 50%
> Enzyme protector 8%
> Non-reactive ingredients 34%

Additional Information for Healthcare Professionals

Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh whole blood samples only. Professionals may use test strips to test capillary and venous whole blood.

Sample Size: 3.0 µL

Reaction Time: 60 seconds

System Measurement Range: 100 mg/dL ~ 400 mg/dL (2.6 mmol/L ~ 10.4 mmol/L)

Hematocrit Range: 20% ~ 60%

Accuracy

Results: Difference Plot Analysis

Model: GD81 / GD82

Number of subjects (N)=320

TABLE 1 Ratio of samples (capillary whole blood) falling within different bias ranges.

Within ± 5 %	Within ± 10 %	Within ± 15 %	Within ± 20 %
80.3%	99.1%	100%	100%

Measurement range: 129 ~ 399 mg/dL

Linear Regression: y = 0.9777x + 1.9945, R² = 0.9724, (n=640)

User Performance

Results: Difference Plot Analysis

Model: GD81 / GD82

Number of subjects (N)=160

TABLE 2 Ratio of samples (capillary whole blood) falling within different bias ranges.

Within ± 5 %	Within ± 10 %	Within ± 15 %	Within ± 20 %
68.4%	98.8%	100%	100%

Measurement range: 129 ~ 399 mg/dL

Linear Regression: y = 0.9753x + 2.0318, R² = 0.9693, (n=320)

Precision

TABLE 3 Precision

GD81 / GD82	Total Cholesterol Concentration (mg/dL)		
	LOW	MID	HIGH
Mean	182.6	241.6	356.6
SD	9.331	11.964	18.594
CV%	5.11%	4.95%	5.21%

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